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Nitrate reduction test

a. Reagents

(1) Nitrate substrate broth : store at 2 to 8 °C.

⎧	peptone	20 g
	(or heart infusion broth for fastidious organisms	25 g)
	potassium nitrate	2 g
	distilled water	1,000 ml

Dispense 2 ml in 13 by 100 mm screw-cap tubes.

Autoclave at 121 °C for 15 min.

(2) 0.8% sulfanilic acid (reagent A)

⎧	sulfanilic acid	0.8 g
	distilled water	70 ml
	glacial acetic acid	30 ml

Mix sulfanilic acid with water ; heat to dissolve.

Cool, and then add acetic acid.

Store at 2 to 8 °C.

Shelf life is 3 months, approximately.

(3) 0.5% N,N-dimethyl- α -naphthyl-amine (reagent B)

⎧	glacial acetic acid	30 ml
	distilled water	70 ml
	MN-dimethyl- α -naphthylamine	0.5 g

Combine acetic acid and water then add α -naphthylamine.

Store at 2 to 8 °C.

Shelf life is 3 months approximately.

(4) Zinc metal dust

b. Procedure

(1) Inoculate the tube of nitrate broth from an isolated colony, a pure subculture, or 1 or 2 drops of an overnight broth culture of the organism.

(2) Incubation

a. Nonfermenting Gram-negative rods, 25 to 30 °C.

b. Other organisms, 35 °C

c. Two to five days of incubation.

d. Incubate *Campylobacter* at 35°C in a microaerobic atmosphere for 72h.

(3) Remove approximately 0.5 ml of broth into a nonsterile 13- by 100-mm tube.

(4) Add 2 or 3 drops of reagent A. Mix well by tapping or shaking tube.

(5) Then, add 2 or 3 drops of reagent B. Mix again.

- (6) Look for a red color within 1 to 2 min.
- (7) If no red color is observed, add a small amount of zinc dust to the nitrate tube.
- (8) Examine for red color within 10 min.

c. Results

Positive results ;

Red color after the addition of reagents: nitrate reduction positive, nitrite reduction negative.

No red color after the addition of reagents plus no red color after the addition of zinc to nitrate broth: nitrate reduction positive, nitrite reduction positive

Negative results ;

In nitrate broth no color development after adding reagents and red color development after adding zinc (zinc catalyzes the change from nitrate to nitrite) : nitrate reduction negative

d. Limitation

Interpretation of color reactions should be made immediately, as color reactions with a positive test may fade rapidly.

(Reference)

CMPH 4th edi. 2016. ASM