RIMD 性状検査法 2022/12/23 版
Biochemical Procedures for the Identification of Bacteria

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Catalase test

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Nitrate Reduction Test

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
【 Reagents 】

Hydrogen peroxide reagent 3%.

Dilute 30% 1:10 in deionized water prior to use.

⇒ 30% reagent is hazardous.

Store at 2 to 8 °C with light shielding.

【 Supplies 】

Glass slide.

Sterile wooden sticks or plastic or platinum loops or wires.

【 Procedure 】

1) Fresh colonies of bacteria growing on agar media, preferably blood agar or Chocolate agar.

For anaerobes, expose colonies to air for 30 min prior to testing.

2) Touch the center of a well isolated colony to a clean glass slide.

⇒ Be sure colony is visible to the naked eye on slide.

⇒ If colony is from blood agar plate, use care not to pick up blood.

3) Place 1 drop of peroxide reagent on slide and observe immediately for effervescence.

⇒ Hold over dark background to well observed bubbles.

4) Discard slide into sharps container.

【 Results 】

Positive results ;

Shows immediate appearance of bubbles or weak reaction has one or two bubbles.

Negative test ;

Shows no bubbles or a few bubbles after 20 sec.

【 Quality control 】

Positive control ;

Bacillus spp. are catalase positive.

Negative control ;

Clostridium spp. are catalase negative.

【 Limitation 】

RBCs contain catalase, do not pick up blood agar with colony.

Selecting colonies with loop will yield false-positive results.

Older cultures may give false-negative results.

Do not mix the reagent and the colony.

Nitrate reduction test

【 Reagents 】

1. Nitrate substrate broth : store at 2 to 8 °C.

[Bacto peptone	20 g
	potassium nitrate	2 g
	distilled water	1,000 ml

Dispense 2 to 4 ml in 13 by 100 mm screw-cap tubes.

Autoclave at 121 °C for 15 min.

2. 0.8% sulfanilic acid (Reagent A)

[sulfanilic acid	0.8 g
	distilled water	70 ml
	glacial acetic acid	30 ml

Mix sulfanilic acid with water ; heat to dissolve.

Cool, and then add acetic acid.

Store at 2 to 8 °C.

Shelf life is 3 months, approximately.

3. 0.5% N,N-dimethyl- α -naphthyl-amine (Reagent B)

[glacial acetic acid	30 ml
	distilled water	70 ml
	MN-dimethyl-1-naphthylamine	0.5 g

Combine acetic acid and water then add α -naphthylamine.

Store at 2 to 8 °C.

Shelf life is 3 months approximately.

4. Zinc metal dust

【 Procedure 】

1) Inoculate the tube of nitrate broth from an isolated colony, a pure subculture, or 1 or 2 drops of an overnight broth culture of the organism.

2) Incubate for 2 hours to 5 days.

Nonfermenting Gram-negative rods incubate at 25 to 30 °C.

Other organisms at 35 to 37 °C.

Incubate *Campylobacter* at 35°C in a microaerobic atmosphere for 3 days.

3) Remove approximately 0.5 ml of broth into a nonsterile 13- by 100-mm tube.

- 4) Add 2 or 3 drops of reagent A. Mix well by tapping or shaking tube.
- 5) Then, add 2 or 3 drops of reagent B. Mix again.
- 6) Look for a red color within 1 to 2 min.
- 7) If no red color is observed, add a small amount of zinc dust to the nitrate tube.
- 8) Examine for red color within 10 min.

【 Results 】

Positive results ;

Red color after the addition of reagents: nitrate reduction positive, nitrite reduction negative.

No red color after the addition of reagents and no red color after the addition of zinc to nitrate broth: nitrate reduction positive, nitrite reduction positive

Negative results ;

In nitrate broth no color development after adding reagents and red color development after adding zinc (zinc catalyses the change from nitrate to nitrite): nitrate reduction negative

【 Limitation 】

Interpretation of color reactions should be made immediately, as color reactions with a positive test may fade rapidly.

【 Reference 】

CMPH 4th edi. 2016. ASM

【 Figure 】